2P *ClearMesh Composite*

IT	ISTRUZIONI PER L'USO
EN	INSTRUCTIONS FOR USE
FR	MODE D'EMPLOI
DE	BEDIENUNGSANLEITUNG
ES	INSTRUCCIONES PARA EL USO
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
PT	INSTRUÇÕES DE UTILIZAÇÃO
PL	INSTRUKCJE OBSŁUGI
RU	ИНСТРУКЦИИ ПО ПРИМЕНЕНИЮ
SR	UPUTSTVA ZA UPOTREBU
HU	HASZNÁLATI UTASÍTÁSOK
TR	KULLANMA TALİMATLARI
AR	إرشادات الاستخدام
CS	NÁVOD K POUŽITÍ
RO	INSTRUCȚIUNI DE UTILIZARE
SK	POKYNÝ K POUŽÍVANIU
BG	ИНСТРУКЦИИ ЗА УПОТРЕБА
ZH	使用說明
NL	GEBRUIKSAANWIJZING
FA	دستور العمل نحوه مصرف

DIPRO MEDICAL DEVICES S.r.l

labeled DIPROMED S.r.l

Via Ciriè 22/a - 10099 San Mauro T.se (TO) - Italy

tel. + 39 011 822 0298 / + 39 011 669 3010

fax. + 39 011 822 3557

www.dipromed.eu - e.mail: info@dipromed.eu

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INSTRUCTIONS FOR USE

COMPOSITE SURGICAL PROSTHESIS IN POLYPROPYLENE MONOFILAMENT, COUPLED TO A POLYPROPYLENE FILM.

Read all instructions carefully. Failure to observe the following instructions may lead to incorrect use of the device and/or injuries to the patient. These instructions are generic instructions for the use of composite surgical prostheses in polypropylene monofilament, coupled to a polypropylene film for intraperitoneal use and must not be intended as a complete reference for the surgical technique to be adopted for the treatment of abdominal wall defects. The device is for the exclusive use of surgeons who have received suitable training for the treatment of abdominal wall defects. The procedure may vary depending on the adopted technique and the anatomical characteristics of the patient.

PATIENT FACTORS

The surgeon must use their experience, good surgical practices to establish whether the 2P composite surgical prosthesis device in polypropylene monofilament, coupled to a polypropylene film, is suited to certain patients. Specific factors concerning patients may compromise healing of the wound with an increased possibility of adverse reactions.

Surgical risks not directly associated with use of the product, its functionality or method of implantation must be understood and evaluated by a qualified surgeon familiar with the anatomy, medical treatments and procedures for the proper treatment of the specific conditions that may be treated using this product. The doctor must provide the patient with accurate and detailed information regarding the possible risks and complications that may arise from the implantation and the adopted surgical procedure.

During the post-operative phase, the surgeon must advise the patient of actions that should be avoided and the possible precautions that should be taken during normal daily activities.

The patient must contact the surgeon immediately in case of fever, reddening and or swelling of the wound, bleeding and any other new symptom that appears following surgery.

DESCRIPTION

Non-absorbable composite prosthesis in polypropylene.

The composite prosthesis is obtained by joining one layer of mesh in polypropylene monofilament and one layer of clear, non-absorbable polypropylene film, for intraperitoneal use.

INDICATIONS

Prosthesis 2P-Cmc - ClearMesh Composite

Composite prosthesis for the treatment of abdominal wall defects, with intraperitoneal positioning. It can be used with either a laparotomic or laparoscopic surgical approach.

Prosthesis 2P-Plurimesh Pcmc - Polyvalent ClearMesh Composite

Cuttable composite prosthesis for the treatment of abdominal wall defects caused by parastomal and/or incisional hernias with intraperitoneal positioning. It can be used with either a laparotomic or laparoscopic surgical approach. The PCMC prosthesis can be adapted to the defect and anatomy of the patient.

Prosthesis 2P-FlaPp Cmc - Free Lateral Polypropylene Prosthesis - ClearMesh Composite

Composite prosthesis for the treatment of large wall defects with loss of peritoneal-fascial substance, with intraperitoneal positioning according to the modified Rives Technique.

Prosthesis 2P-UJCMC – Umbilical ClearMesh Composite

Composite prosthesis for the treatment of abdominal wall defects, in particular umbilical and epigastric hernias, small incisional trocar-site hernias, with intraperitoneal or preperitoneal positioning. It can be used with either a laparotomic or laparoscopic surgical approach.

In line with good surgical practice, the intended use of the wings is merely for the final suture of the prosthesis along the edges of the defect. The excess part is to be removed. The positioning wire is to be used only for traction and the correct positioning of the prosthesis. Once positioning is complete, the wire must be removed.

PROCEDURE FOR USE

The device is for the exclusive use of specialised medical personnel, trained in its use.

The product is a non-absorbable implantable medical device and can only be used in operating theatres.

The dimensions of the prosthesis must be determined in accordance with good surgical practices and depending on the size of the defect, any possible narrowing of the prosthesis, the condition and pathology of the host tissue.

The prosthesis must be affixed as per the procedures required for the fixing system and according to best surgical practice.

In the event the device needs to be removed, the surgeon must follow good clinical and surgical practices and ensure relevant procedures are followed in order to guarantee the patient's safety.

USE AND POSITIONING

To guarantee its functionality it is essential that the product be correctly oriented.

The surface in polypropylene film is made for the peritoneal cavity; the film side of the prosthesis can come into contact with the intestine or other visceral structures, taking care not to damage the surface in polypropylene film.

The polypropylene monofilament mesh side has the same colonisation characteristics as single-layer polypropylene meshes. Therefore, this surface must be facing the part on which tissue growth is sought to be promoted. The surface of the prosthesis in polypropylene mesh must never be positioned in contact with the intestine or visceral structures.

Only the PCMC prosthesis can be adapted to size. Use a sharpened tool to cut the prosthesis. The distribution of the seams on the prosthesis allows the prosthesis to be cut without separating the two layers.

CONTRAINDICATIONS

The polypropylene monofilament mesh side must never be positioned in contact with the intestine or visceral structures, as this may cause serious adhesions.

Specific factors concerning patients may compromise healing of the wound with an increased possibility of adverse reactions.

Do not use the prosthesis in growing patients, in patients who are pregnant or who are planning a pregnancy.

Do not use on infected and contaminated tissues.

WARNINGS

- The device is for the exclusive use of specialised medical personnel, trained in its use;
- the device is supplied sterile. Check the packaging and make sure the product is intact before use;
- the doctor and his/her collaborators must visually inspect the packaged device to check that the sterile barrier has not been damaged;
- do not use if the package is open or damaged, insofar as the product's sterility may have been compromised;
- do not use the device after the expiry date;
- the device is disposable and is intended to be used once only, for one patient only, insofar as it cannot be re-sterilised and is non-reusable. Reprocessing may compromise the stability of the chemical-physical-mechanical and biological properties of the device;
- guarantee correct orientation and positioning: position the film side of the prosthesis in the peritoneal cavity, orienting it toward the intestine or sensitive organs, taking care not to damage the surface in polypropylene film. Do not place the mesh surface in contact with visceral structures, in particular with intestinal loops;
- in case of infection, the prosthesis may need to be removed;
- complications relative to use of the prosthesis may include intestinal obstruction, adhesions, extrusion, migrations, the formation of fistulas, relapses, seromas, haematomas and inflammations, infections;
- the user must evaluate the correct size of the prosthesis in relation to the defect to prevent any relapses;
- use of any prosthesis in a contaminated context may lead to the formation of fistulas and/or expulsion of the prosthesis;
- during the introduction and positioning of the device, manoeuvres must be carried out taking care not to damage the device;
- the patient must contact the surgeon immediately in case of fever, reddening and or swelling of the wound, bleeding and any other new symptom that appears following surgery.

PRECAUTIONS

- Read all the instructions carefully before use;
- open the package in the operating theatre and in accordance with good practices;
- once the sterile packaging of the prosthesis has been opened, it cannot be re-sterilised or reused;
- the prosthesis must be stored in a clean, dry place far from sources of light and heat;
- handle the product with sterile gloves and atraumatic instrumentation;
- medical staff must act in accordance with good clinical practices and instruct the patient during the post-operative phase regarding actions to be avoided and precautions to be taken during normal daily activities;
- after the surgery and for a minimum of 30 days, the patient must avoid lifting weights, strain and sports activities involving the abdominal wall. Normal activities may be resumed after a check-up and authorisation by the doctor that performed the surgery.

ADVERSE REACTIONS

Possible adverse reactions include blockages, adhesions, extrusion, migration, formation of fistulas, relapses, seromas, haematomas and inflammations, infections and relapse of the hernial pathology. Other known adverse reactions relative to surgical procedures include pain, bleeding, failure of the procedure;

As with all implantations, local inflammation may occur in the wound site and/or a foreign body reaction;

As with all foreign bodies, the polypropylene prosthesis may cause the clinical worsening of an existing infection; in these circumstances, it may be necessary to remove the device;

STERILISATION

All composite prostheses are disposable and sterilised by ethylene oxide (EO).
Do not use the device if the package is damaged.

TRACEABILITY

Identification labels are included in each package for product traceability, as provided for by legislation in force.

DISPOSAL

After use, eliminate the product and its relative packaging in compliance with hospital, administrative and/or local regulations in force

EXPIRY DATE

The product has a 5-year expiry date.

EXPECTED LIFE

For functional use the implanted device has a long-term life.

PACKAGING

The commercial box constitutes a unit of sale; the product inside is not sold separately.
Any infusions in the packaging do not constitute mobile particulate matter

THE PRODUCT IS LATEX-FREE, PHTHALATE-FREE, NON-PYROGENIC

In case of disputes regarding interpretation and/or the deviation of content, the Italian text shall prevail.

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